

Centeze™



Centeze™ Centesis Catheter

Galt Medical Corp. is pleased to introduce the Centeze™ over the needle catheter system for centesis procedures. Featuring large skived holes, a smooth transition, and secure slip fit needle design, the Galt Centeze was designed to enhance ease in every centesis procedure.

- 3 large skived holes to aid in maximum drainage and reduction of clogging
- Tapered tip, smooth transition, and siliconized catheter to aid in insertion
- Echogenic needle with bevel indicator to aid in visibility
- Clear needle hub for flashback confirmation




For more information about Centeze™ Centesis Catheter or to set up a meeting with a Galt Sales Consultant call (800) 639-2800.

Centeze™



Clear needle hub



Large skived holes



Smooth Transition

Centeze™ Centesis Catheter

Catheter	Catheter Size	Catheter Length	Needle Gauge	Needle Length	Needle Tip	Lock	Box Quantity
DRC-002-01	4Fr	7cm	21ga	9cm	Echogenic	Slip fit	10
DRC-002-04	4Fr	10cm	21ga	12cm	Echogenic	Slip fit	10
DRC-002-02	5Fr	7cm	19ga	9cm	Echogenic	Slip fit	10
DRC-002-05	5Fr	10cm	19ga	12cm	Echogenic	Slip fit	10
DRC-002-08	5Fr	15cm	19ga	17cm	Echogenic	Slip fit	10

