

## Fine Needle Aspiration (FNA) Biopsy Needles

Our traditional FNA biopsy needles harvest aspirate for cytological evaluations or specimens for histological evaluations.

### Key Product Attributes

- Each needle has numerically ordered centimeter markings on the outer cannula to provide reference for depth placement
- Needles feature an adjustable needle stop which allows the user to restrict forward movement, localizing the needle tip to the biopsy site
- Color coded stylet hubs indicate gauge size of needles
- Available in a variety of gauge sizes and centimeter lengths
- Available with an echogenic tip on the outer cannula to promote accurate placement under ultrasound guidance

### Chiba Biopsy Needles (Box of 10)

CATALOG NO.	DESCRIPTION
MCN1803	18G x 9cm
MCN1803US	18G x 9cm with echogenic tip
MCN1806	18G x 15cm
MCN1806US	18G x 15cm with echogenic tip
MCN1808	18G x 20cm
MCN1808US	18G x 20cm with echogenic tip
MCN2003	20G x 9cm
MCN2003US	20G x 9cm with echogenic tip
MCN2006	20G x 15cm
MCN2006US	20G x 15cm with echogenic tip
MCN2008	20G x 20cm
MCN2008US	20G x 20cm with echogenic tip
MCN2203	22G x 9cm
MCN2203US	22G x 9cm with echogenic tip
MCN2206	22G x 15cm
MCN2206US	22G x 15cm with echogenic tip
MCN2208	22G x 20cm
MCN2208US	22G x 20cm with echogenic tip

### Franseen Biopsy Needles (Box of 10)

CATALOG NO.	DESCRIPTION
MFN1803	18G x 9cm
MFN1803US	18G x 9cm with echogenic tip
MFN1806	18G x 15cm
MFN1806US	18G x 15cm with echogenic tip
MFN1808	18G x 20cm
MFN1808US	18G x 20cm with echogenic tip
MFN2003	20G x 9cm
MFN2003US	20G x 9cm with echogenic tip
MFN2006	20G x 15cm
MFN2006US	20G x 15cm with echogenic tip
MFN2008	20G x 20cm
MFN2008US	20G x 20cm with echogenic tip
MFN2203	22G x 9cm
MFN2203US	22G x 9cm with echogenic tip
MFN2206	22G x 15cm
MFN2206US	22G x 15cm with echogenic tip
MFN2208	22G x 20cm
MFN2208US	22G x 20cm with echogenic tip



### Spinal Biopsy Needles (Box of 10)

CATALOG NO.	DESCRIPTION
MSPN1803	18G x 9cm
MSPN1806	18G x 15cm
MSPN2003	20G x 9cm
MSPN2003US	20G x 9cm with echogenic tip
MSPN2006	20G x 15cm
MSPN2006US	20G x 15cm with echogenic tip
MSPN2203	22G x 9cm
MSPN2203US	22G x 9cm with echogenic tip
MSPN2206	22G x 15cm
MSPN2206US	22G x 15cm with echogenic tip
MSPN2208	22G x 20cm
MSPN2502	25G x 4cm
SPN2302	23G x 5cm
SPN2303	23G x 9cm
SPN2502	25G x 5cm
SPN2503	25G x 9cm

### Westcott Biopsy Needles (Box of 10)

CATALOG NO.	DESCRIPTION
MWN2003	20G x 9cm
MWN2003US	20G x 9cm with echogenic tip
MWN2006	20G x 15cm
MWN2006US	20G x 15cm with echogenic tip
MWN2008	20G x 20cm
MWN2008US	20G x 20cm with echogenic tip
MWN2203	22G x 9cm
MWN2203US	22G x 9cm with echogenic tip
MWN2206	22G x 15cm
MWN2206US	22G x 15cm with echogenic tip
MWN2208	22G x 20cm
MWN2208US	22G x 20cm with echogenic tip

# Soft Tissue Biopsy

## Fine Needle Aspiration (FNA) Biopsy Needles (cont.)

Trocar Biopsy Needles (Box of 10)	
CATALOG NO.	DESCRIPTION
MTN1806	18G x 15cm
MTN1808	18G x 20cm

Greene Biopsy Needles (Box of 10)	
CATALOG NO.	DESCRIPTION
MGN2206X	22G x 15cm includes 19G x 10cm introducer with echogenic tip

10cc Locking Syringe (Box of 10)	
CATALOG NO.	DESCRIPTION
MLS1001	10cc

## Fine Needle Aspiration Trays

	16G x 10cm Techna-Cut™ Biopsy Needle (Box of 10)	18G x 10cm Techna-Cut™ Biopsy Needle (Box of 10)	16G x 7cm Techna-Cut™ Biopsy Needle (Box of 10)	18G x 7cm Techna-Cut™ Biopsy Needle (Box of 10)
Trays include	ASK1016A	ASK1017A	ASK7016A	ASK7017A
25G x 5/8" (1.6cm) Needle	1	1	1	1
21G x 1 1/2" (3.8cm) Needle	1	1	1	1
20G x 1 1/2" (3.8cm) Needle	1	1	1	1
5mL Luer-lock Syringe	1	1	1	1
Thumb Scalpel	1	1	1	1
Amber Specimen Vials	1	1	1	1
3"x 3"Gauze Pads	3	3	3	3
10% Iodophor PVP Swab Stick	3	3	3	3
1% Lidocaine HCL Ampule (5mL)	1	1	1	1
10mL Sodium Chloride	1	1	1	1
Label Sheet	1	1	1	1
Fenestrated Drape	1	1	1	1
Towel	1	1	1	1
Bandage	1	1	1	1
Hospital Wrap	1	1	1	1