



Access with Ease

V•Stick Vascular Access Sets are designed to facilitate a smooth and easy introduction into the vessel for guidewires up to .038" and minimize trauma to the surrounding tissue.


Access Needle



Coaxial Introducer



Nitinol or Stainless Steel Guidewire





Ordering Information

V•Stick Vascular Access Sets w/Standard Co-axial Introducer (Box of 10)							
CATALOG #	SHEATH SIZE	NEEDLE GAUGE x LENGTH	ECHOGENIC	STIFFNESS	DILATOR LOCKING NUT COLOR	GUIDEWIRE DIAMETER x LENGTH	MATERIAL
390204400E	4F	21G x 3.8cm	No	Standard	Red	.018" x 40cm	NT
391504000E	4F	21G x 7cm	Yes	Standard	Red	.018" x 40cm	SS
391504300E	4F	21G x 7cm	Yes	Standard	Red	.018" x 40cm	NT
391504400E	4F	21G x 7cm	No	Standard	Red	.018" x 40cm	NT
391504500E	4F	21G x 7cm	No	Standard	Red	.018" x 40cm	SS
391505000E	5F	21G x 7cm	Yes	Standard	Red	.018" x 40cm	SS
391505300E	5F	21G x 7cm	Yes	Standard	Red	.018" x 40cm	NT
391505400E	5F	21G x 7cm	No	Standard	Red	.018" x 40cm	NT
391505500E	5F	21G x 7cm	No	Standard	Red	.018" x 40cm	SS

V•Stick Vascular Access Sets w/Stiff Co-axial Introducer (Box of 10)							
CATALOG #	SHEATH SIZE	NEEDLE GAUGE x LENGTH	ECHOGENIC	STIFFNESS	DILATOR LOCKING NUT COLOR	GUIDEWIRE DIAMETER x LENGTH	MATERIAL
391604000E	4F	21G x 7cm	Yes	Stiff	White	.018" x 40cm	SS
391604300E	4F	21G x 7cm	Yes	Stiff	White	.018" x 40cm	NT
391604400E	4F	21G x 7cm	No	Stiff	White	.018" x 40cm	NT
391604500E	4F	21G x 7cm	No	Stiff	White	.018" x 40cm	SS
391605000E	5F	21G x 7cm	Yes	Stiff	White	.018" x 40cm	SS
391605300E	5F	21G x 7cm	Yes	Stiff	White	.018" x 40cm	NT
391605400E	5F	21G x 7cm	No	Stiff	White	.018" x 40cm	NT
391605500E	5F	21G x 7cm	No	Stiff	White	.018" x 40cm	SS

SS = Stainless Steel

NT = Nitinol