

MICRO-INTRODUCER KITS & STIFFEN MICRO-INTRODUCER KITS




Galt's Micro-Introducer Kits provide the high-quality and precision design expected in a fine medical device. Each feature was designed to aid in the access process, allowing for consistent performance every time.


Galt's line of Stiffen Micro-Introducers sets a new standard for high-quality and precision design in micro-introducers. Kits combine our precise dilator and sheath transitions with a reinforced assembly to facilitate smooth entry into difficult vessels.

4 and 5 French sizes available
Variety of wire configurations

Standard or Echogenic needles
All-inclusive packaging



Smooth transition



Kink-resistant core-wire available



Needle bevel indicator



For more information about the Micro-Introducers or to set up a meeting with a Galt Sales Consultant call (800) 639-2800.

Micro-Introducer Kits

Catalog Number	Sheath Size	Sheath Length	Guidewire	Guidewire Type	Needle Length	Needle Type
KIT-002-30	4.0F	10cm	40cm	SS/SS	21ga x 7cm	Echogenic
KIT-002-31	5.0F	10cm	40cm	SS/SS	21ga x 7cm	Echogenic
KIT-002-36	4.0F	10cm	60cm	SS/SS	21ga x 7cm	Echogenic
KIT-002-37	5.0F	10cm	60cm	SS/SS	21ga x 7cm	Echogenic
KIT-002-28	4.0F	10cm	40cm	SS/Galt Gold®	21ga x 7cm	Standard
KIT-002-29	5.0F	10cm	40cm	SS/Galt Gold®	21ga x 7cm	Standard
KIT-002-64	5.0F	10cm	60cm	Nitinol/SS	21ga x 4cm	Smooth, non-echogenic
KIT-002-32	4.0F	10cm	40cm	Nitinol/SS	21ga x 7cm	Echogenic
KIT-002-33	5.0F	10cm	40cm	Nitinol/SS	21ga x 7cm	Echogenic
KIT-002-38	4.0F	10cm	60cm	Nitinol/SS	21ga x 7cm	Echogenic
KIT-002-39	5.0F	10cm	60cm	Nitinol/SS	21ga x 7cm	Echogenic
KIT-002-12	4.0F	10cm	40cm	Nitinol/Galt Gold®	21ga x 7cm	
KIT-002-13	5.0F	10cm	40cm	Nitinol/Galt Gold®	21ga x 7cm	
KIT-002-34	4.0F	10cm	40cm	Nitinol/Galt Gold®	21ga x 7cm	Echogenic
KIT-002-35	5.0F	10cm	40cm	Nitinol/Galt Gold®	21ga x 7cm	Echogenic
KIT-002-62	4.0F	10cm	40cm	Nitinol/Galt Gold®	21ga x 4cm	
KIT-002-63	5.0F	10cm	40cm	Nitinol/Galt Gold®	21ga x 4cm	
KIT-002-43	5.0F	10cm	40cm	Nitinol/SS	21ga x 7cm	Echogenic

Stiffen Micro-Introducer Kits

Catalog Number	Sheath Size	Sheath Length	Guidewire	Guidewire Type	Needle Length	Needle Type
KIT-011-02	4.0F	10cm	40cm	SS/SS	21ga x 7cm	Echogenic
KIT-011-03	5.0F	10cm	40cm	SS/SS	21ga x 7cm	Echogenic
KIT-011-40	4.0F	10cm	40cm	Nitinol/SS	21ga x 7cm	Standard
KIT-011-41	5.0F	10cm	40cm	Nitinol/SS	21ga x 7cm	Standard
KIT-011-42	4.0F	10cm	40cm	Nitinol/SS	21ga x 7cm	Echogenic
KIT-011-43	5.0F	10cm	40cm	Nitinol/SS	21ga x 7cm	Echogenic
KIT-011-62	4.0F	10cm	40cm	Nitinol/Galt Gold®	21ga x 7cm	Echogenic
KIT-011-63	5.0F	10cm	40cm	Nitinol/Galt Gold®	21ga x 7cm	Echogenic

Micro-Introducer Kits with Side Holes

Catalog Number	Sheath Size	Sheath Length	Guidewire	Guidewire Type	Needle Length	Needle Type
KIT-033-04	4.0F	10cm	40cm	Nitinol/SS	21ga x 7cm	Echogenic

Bariatric Micro-Introducer Kits

Catalog Number	Sheath Size	Sheath Length	Guidewire	Guidewire Type	Needle Length	Needle Type
KIT-011-99	4.0F	15cm	60cm	Nitinol/Galt Gold®	21ga x 12cm	Echogenic



NOTE: Order Quantity – 10 kits/box