

## Non-Positionable Breast Localization Needles (BLN)

### Hawkins™ III BLN

The Hawkins™ III BLN offers a barb retention design to provide strong anchoring strength in the breast.

#### Key Product Attributes

- Available in two wire choices: FlexStrand™, a 19-strand cable providing maximum flexibility with minimal risk of transection, or a traditional hardwire that minimizes flex during placement
- Skin retention clip prevents forward migration of wire during patient transport
- Wire is marked in two places: I indicates wire at distal tip of needle; II indicates barb deployment
- Standard marking cannula on wire, beginning 2.5 cm from distal tip, serves as a landmark for the surgeon
- End hole design allows for dye injection or fluid aspiration
- Centimeter markings provide reference for depth placement
- Available with echogenic tip

#### Hawkins™ III Breast Localization Needles w/Echogenic Tip and Hardwire (Box of 10)

CATALOG NO.	DESCRIPTION
244030	20G x 3 cm, hardwire with echogenic tip
244050	20G x 5 cm, hardwire with echogenic tip
244075	20G x 7.5 cm, hardwire with echogenic tip
244100	20G x 10 cm, hardwire with echogenic tip
244125	20G x 12.5 cm, hardwire with echogenic tip

#### Hawkins™ III Breast Localization Needles w/FlexStrand™ (Box of 10)

CATALOG NO.	DESCRIPTION
253030	20G x 3 cm
253050	20G x 5 cm
253075	20G x 7.5 cm
253100	20G x 10 cm
253125	20G x 12.5 cm

#### Hawkins™ III Breast Localization Needles w/Hardwire (Box of 10)

CATALOG NO.	DESCRIPTION
243030	20G x 3 cm
243050	20G x 5 cm
243075	20G x 7.5 cm
243100	20G x 10 cm
243125	20G x 12.5 cm

