

Guidewire Introducer Needles

Key Product Attributes

- Designed for anterior, single wall, arterial percutaneous puncture
- Tapered luer hub allows for easy guidewire insertion
- Available with a clear window wall hub
- Also available with winged base plate

Guidewire Introducer Needle (Box of 25)

CATALOG NO.	DESCRIPTION
GW1802	18G x 7cm, accepts .038" guidewire
GW1902	19G x 7cm, accepts .035" guidewire
GW1802WW	18G x 7cm, accepts .038" guidewire, window wall hub

Guidewire Introducer Needle/Winged Base Plate (Box of 25)

CATALOG NO.	DESCRIPTION
BGW1802	18G x 7cm, accepts .038" guidewire
BGW1902	19G x 7cm, accepts .035" guidewire

Seldinger Needles

- Designed for easy, atraumatic, percutaneous arterial puncture
- Color coded stylet

Seldinger Needle (Box of 25)

CATALOG NO.	DESCRIPTION
SN1802	18G x 7cm, accepts .038" guidewire
SN1902	19G x 7cm, accepts .035" guidewire

Potts/Cournand Needle

- Designed for carotid angiography, direct arterial pressure monitoring, blood sampling and percutaneous catheterization

Potts/Cournand Needle (Box of 25)

CATALOG NO.	DESCRIPTION
MPCN1802	18G x 5.7cm, accepts .038" guidewire

Vessel Dilators

Feathered Tip

CATALOG NO.	SIZE	LENGTH	FOR GUIDEWIRE SIZES	UNITS
392354	4F	19.8cm	Accepts .018" - .021"	10
392355	5F	19.8cm	Accepts .018" - .021"	10
392356	6F	19.8cm	Accepts .018" - .021"	10

Standard Taper

CATALOG NO.	SIZE	LENGTH	FOR GUIDEWIRE SIZES	UNITS
392299	4F	13.2cm	Accepts .032" - .035"	10
392310	5F	19.8cm	Accepts .035" - .038"	10
392305	6F	19.8cm	Accepts .035" - .038"	10
392306	7F	19.8cm	Accepts .035" - .038"	10
392307	8F	19.8cm	Accepts .035" - .038"	10
392309	9F	19.8cm	Accepts .035" - .038"	10

